

The Multi-Lift can be used indoors in bedrooms, bathrooms or any room!
The Lift is powed by either a battery pack or a conveinent transformer.

A handy Magnetic control box attached to a spiral cord, allows fleible transfers and is easily stored on any metal surface. Sidways movents are done maunually.

For use in the Home, simply install a wall or floor mount where the lift is to be used. The lift can be easily

adjusted for different heights by adjusting the braket on the mount.

Many lift extensions and adapters are available to cover the widest possible area or make lifting easy, even when the area is tight.

Available in either the standard 275lb lifting capacity or in a heavy duty version that will carry 400lbs, the Multi-Lift is perfect for most any patient.

The Multi-Lift is compact as well, weighing in at just 23lbs, the Multi-lift is easily lifted off its mounting bracket when not in use. Its small size, only 22.5 inches long and 4 inches wide, the lift is easily stored.















- Powered by a conveinent transformer at each wallmount station
- Used for lifting the person from the wheelchair to the bed, toilet, bath or ANY area
- Same lift is used for ALL mounts!
- Numerous extensions and adapters available



1-800-849-2143 www.gccessunlimited.com